




VDAC1 Antibody

Product Code	CSB-PA005626
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21796
Immunogen	Synthesized peptide derived from the N-terminal region of Human VDAC1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VDAC1; VDAC; Voltage-dependent anion-selective channel protein 1; VDAC-1; hVDAC1; Outer mitochondrial membrane protein porin 1; Plasmalemmal porin; Porin 31HL; Porin 31HM
Immunogen Species	Homo sapiens (Human)
Target Names	VDAC1

Image



Western Blot analysis of K562 cells using VDAC1 Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.