




GP1BA Antibody

Product Code	CSB-PA005629
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07359
Immunogen	Synthesized peptide derived from the Internal region of Human CD42b.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GP1BA; Platelet glycoprotein Ib alpha chain; GP-Ib alpha; GPIb-alpha; GPIbA; Glycoprotein Ibalpha; Antigen CD42b-alpha; CD42b
Immunogen Species	Homo sapiens (Human)
Target Names	GP1BA

Image



Western Blot analysis of HeLa cells using CD42b Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.