




MRC2 Antibody

Product Code	CSB-PA005646
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBG0
Immunogen	Synthesized peptide derived from the Internal region of Human Endo180.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRC2; CLEC13E; ENDO180; KIAA0709; UPARAP; C-type mannose receptor 2; C-type lectin domain family 13 member E; Endocytic receptor 180; Macrophage mannose receptor 2
Immunogen Species	Homo sapiens (Human)
Target Names	MRC2

Image



Western Blot analysis of HeLa cells using Endo180 Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.