




CD163L1 Antibody

Product Code	CSB-PA005658
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR16
Immunogen	Synthesized peptide derived from the Internal region of Human CD163b.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD163L1; CD163B; M160; Scavenger receptor cysteine-rich type 1 protein M160; CD163 antigen-like 1; CD163b
Immunogen Species	Homo sapiens (Human)
Target Names	CD163L1

Image



Western Blot analysis of HeLa cells using CD163b Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.