




INHBA Antibody

Product Code	CSB-PA005668
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08476
Immunogen	Synthesized peptide derived from the C-terminal region of Human Inhibin β -A.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	INHBA; Inhibin beta A chain; Activin beta-A chain; Erythroid differentiation protein; EDF
Immunogen Species	Homo sapiens (Human)
Target Names	INHBA

Image



Western Blot analysis of MCF-7, 293 cells using Inhibin β -A Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.