




THPO Antibody

Product Code	CSB-PA005677
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40225
Immunogen	Synthesized peptide derived from the Internal region of Human TPO.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	THPO; MGDF; Thrombopoietin; C-mpl ligand; ML; Megakaryocyte colony-stimulating factor; Megakaryocyte growth and development factor; MGDF; Myeloproliferative leukemia virus oncogene ligand
Immunogen Species	Homo sapiens (Human)
Target Names	THPO

Image



Western Blot analysis of MCF-7 cells using TPO Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.