



KLF1 Antibody

Product Code	CSB-PA005718
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13351
Immunogen	Synthesized peptide derived from Human EKLF around the non-acetylation site of K274.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA; ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KLF1; EKLF; Krueppel-like factor 1; Erythroid krueppel-like transcription factor; EKLF
Immunogen Species	Homo sapiens (Human)
Target Names	KLF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.