




MYB Antibody

Product Code	CSB-PA005729
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10242
Immunogen	Synthesized peptide derived from Human c-Myb around the non-acetylation site of K480.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYB

Image



Western Blot analysis of MCF7 cells using c-Myb Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.