




KAT2B Antibody

Product Code	CSB-PA005734
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92831
Immunogen	Synthesized peptide derived from Human PCAF around the non-acetylation site of K428.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KAT2B; PCAF; Histone acetyltransferase KAT2B; Histone acetyltransferase PCAF; Histone acetylase PCAF; Lysine acetyltransferase 2B; P300/CBP-associated factor; P/CAF
Immunogen Species	Homo sapiens (Human)
Target Names	KAT2B

Image



Western Blot analysis of PC12 cells using PCAF Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.