




FOXP3 Antibody

Product Code	CSB-PA005784
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZS1
Immunogen	Synthesized peptide derived from the C-terminal region of Human FOXP3.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXP3; IPEX; JM2; Forkhead box protein P3; Scurfin
Immunogen Species	Homo sapiens (Human)
Target Names	FOXP3

Image



Western Blot analysis of PC12 cells using FOXP3 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.