




FMO3 Antibody

Product Code	CSB-PA005810
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31513
Immunogen	Synthesized peptide derived from the Internal region of Human FMO3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FMO3; Dimethylaniline monooxygenase [N-oxide-forming] 3; Dimethylaniline oxidase 3; FMO II; FMO form 2; Hepatic flavin-containing monooxygenase 3; FMO 3; Trimethylamine monooxygenase
Immunogen Species	Homo sapiens (Human)
Target Names	FMO3

Image



Western Blot analysis of HeLa cells using FMO3 Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.