




CD74 Antibody

Product Code	CSB-PA005842
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04233
Immunogen	Synthesized peptide derived from the Internal region of Human CD74.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD74; DHLAG; HLA class II histocompatibility antigen gamma chain; HLA-DR antigens-associated invariant chain; Ia antigen-associated invariant chain; Ii; p33; CD74
Immunogen Species	Homo sapiens (Human)
Target Names	CD74

Image



Western Blot analysis of PC12 cells using CD74 Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.