




LILRA5 Antibody

Product Code	CSB-PA005843
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NI73
Immunogen	Synthesized peptide derived from the Internal region of Human CD85f.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LILRA5; ILT11; LILRB7; LIR9; Leukocyte immunoglobulin-like receptor subfamily A member 5; CD85 antigen-like family member F; Immunoglobulin-like transcript 11; ILT-11; Leukocyte immunoglobulin-like receptor 9; LIR-9; CD85f
Immunogen Species	Homo sapiens (Human)
Target Names	LILRA5

Image



Western Blot analysis of HeLa cells using CD85f Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.