




IL18R1 Antibody

Product Code	CSB-PA005855
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13478
Immunogen	Synthesized peptide derived from the Internal region of Human IL-18R.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL18R1; IL1RRP; Interleukin-18 receptor 1; IL-18R-1; IL-18R1; CD218 antigen-like family member A; CDw218a; IL1 receptor-related protein; IL-1Rrp; IL1R-rp; CD218a
Immunogen Species	Homo sapiens (Human)
Target Names	IL18R1

Image



Western Blot analysis of 293 cells using IL-18R Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.