




CD200 Antibody

Product Code	CSB-PA005858
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41217
Immunogen	Synthesized peptide derived from the Internal region of Human CD200.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD200; MOX1; MOX2; My033; OX-2 membrane glycoprotein; CD200
Immunogen Species	Homo sapiens (Human)
Target Names	CD200

Image



Western Blot analysis of PC12 cells using CD200 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.