




POLR3G Antibody

Product Code	CSB-PA005904
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15318
Immunogen	Synthesized peptide derived from the Internal region of Human POLR3G.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	POLR3G; DNA-directed RNA polymerase III subunit RPC7; RNA polymerase III subunit C7; DNA-directed RNA polymerase III subunit G; RNA polymerase III 32 kDa subunit; RPC32
Immunogen Species	Homo sapiens (Human)
Target Names	POLR3G

Image



Western Blot analysis of PC12 cells using POLR3G Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.