




SAMHD1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA005912 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y3Z3 |
| Immunogen | Synthesized peptide derived from the Internal region of Human SAMHD1. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | SAMHD1; MOP5; SAM domain and HD domain-containing protein 1; Dendritic cell-derived IFNG-induced protein; DCIP; Monocyte protein 5; MOP-5 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SAMHD1 |

Image



Western Blot analysis of PC12 cells using SAMHD1 Polyclonal Antibody.