




HLA-DQB1/HLA-DQB2 Antibody

Product Code	CSB-PA005915
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01920/P05538
Immunogen	Synthesized peptide derived from the Internal region of Human HLA-DQB1/2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HLA-DQB1; HLA-DQB; HLA class II histocompatibility antigen, DQ beta 1 chain; MHC class II antigen DQB1; HLA-DQB2; HLA-DXB; HLA class II histocompatibility antigen, DQ beta 2 chain
Immunogen Species	Homo sapiens (Human)
Target Names	HLA-DQB1/HLA-DQB2

Image



Western Blot analysis of PC12 cells using HLA-DQB1/2 Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.