




EEF1A1/EEF1A2/EEF1A1P5 Antibody

Product Code	CSB-PA005918
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68104/Q05639/Q5VTE0
Immunogen	Synthesized peptide derived from the N-terminal region of Human EF-1 α 1/2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EEF1A1; EEF1A; EF1A; LENG7; Elongation factor 1-alpha 1; EF-1-alpha-1; Elongation factor Tu; EF-Tu; Eukaryotic elongation factor 1 A-1; eEF1A-1; Leukocyte receptor cluster member 7; EEF1A2; EEF1AL; STN; Elongation factor 1-alpha 2; EF-1-alpha-2; Eukaryoti
Immunogen Species	Homo sapiens (Human)
Target Names	EEF1A1/EEF1A2/EEF1A1P5

Image



Western Blot analysis of HepG2 cells using EF-1 α 1/2 Polyclonal Antibody