




SRY Antibody

Product Code	CSB-PA005928
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05066
Immunogen	Synthesized peptide derived from the Internal region of Human SRY.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SRY; TDF; Sex-determining region Y protein; Testis-determining factor
Immunogen Species	Homo sapiens (Human)
Target Names	SRY

Image



Western Blot analysis of HeLa cells using SRY Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.