





CD47 Antibody

Product Code	CSB-PA005993
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08722
Immunogen	Synthesized peptide derived from the Internal region of Human CD47.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD47; MER6; Leukocyte surface antigen CD47; Antigenic surface determinant protein OA3; Integrin-associated protein; IAP; Protein MER6; CD47
Immunogen Species	Homo sapiens (Human)
Target Names	CD47

Image



Western Blot analysis of 3T3 cells using CD47 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using CD47 Polyclonal Antibody