




LILRB5 Antibody

Product Code	CSB-PA005995
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75023
Immunogen	Synthesized peptide derived from the N-terminal region of Human CD85c.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LILRB5; LIR8; Leukocyte immunoglobulin-like receptor subfamily B member 5; CD85 antigen-like family member C; Leukocyte immunoglobulin-like receptor 8; LIR-8; CD85c
Immunogen Species	Homo sapiens (Human)
Target Names	LILRB5

Image



Western Blot analysis of HepG2 cells using CD85c Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.