




TNFSF10 Antibody

Product Code	CSB-PA005996
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50591
Immunogen	Synthesized peptide derived from the Internal region of Human TRAIL.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFSF10; APO2L; TRAIL; Tumor necrosis factor ligand superfamily member 10; Apo-2 ligand; Apo-2L; TNF-related apoptosis-inducing ligand; Protein TRAIL; CD253
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF10

Image



Western Blot analysis of KB cells using TRAIL Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.