




TERT Antibody

Product Code	CSB-PA006014
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14746
Immunogen	Synthesized peptide derived from the C-terminal region of Human TERT.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TERT; EST2; TCS1; TRT; Telomerase reverse transcriptase; HEST2; Telomerase catalytic subunit; Telomerase-associated protein 2; TP2
Immunogen Species	Homo sapiens (Human)
Target Names	TERT

Image



Western Blot analysis of HeLa HepG2 cells using TERT Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.