




FLT3LG Antibody

Product Code	CSB-PA006021
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49771
Immunogen	Synthesized peptide derived from the C-terminal region of Human Flt3-L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FLT3LG; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; Flt3L; SL cytokine
Immunogen Species	Homo sapiens (Human)
Target Names	FLT3LG

Image



Western Blot analysis of HeLa cells using Flt3-L Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.