




SPARC Antibody

Product Code	CSB-PA006034
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09486
Immunogen	Synthesized peptide derived from the Internal region of Human SPARC.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SPARC; ON; SPARC; Basement-membrane protein 40; BM-40; Osteonectin; ON; Secreted protein acidic and rich in cysteine
Immunogen Species	Homo sapiens (Human)
Target Names	SPARC

Image



Western Blot analysis of HeLa HepG2 cells using SPARC Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.