




CD63 Antibody

Product Code	CSB-PA006039
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08962
Immunogen	Synthesized peptide derived from the Internal region of Human CD63.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD63; MLA1; TSPAN30; CD63 antigen; Granulophysin; Lysosomal-associated membrane protein 3; LAMP-3; Melanoma-associated antigen ME491; OMA81H; Ocular melanoma-associated antigen; Tetraspanin-30; Tspan-30; CD63
Immunogen Species	Homo sapiens (Human)
Target Names	CD63

Image



Western Blot analysis of SKOV3 cells using CD63 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.