




CSF3R Antibody

Product Code	CSB-PA006057
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99062
Immunogen	Synthesized peptide derived from the Internal region of Human G-CSFR.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSF3R; GCSFR; Granulocyte colony-stimulating factor receptor; G-CSF receptor; G-CSF-R; CD114
Immunogen Species	Homo sapiens (Human)
Target Names	CSF3R

Image



Western Blot analysis of SKOV3 cells using G-CSFR Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.