




ZNF384 Antibody

Product Code	CSB-PA006074
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TF68
Immunogen	Synthesized peptide derived from the N-terminal region of Human ZNF384.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF384; CAGH1; CIZ; NMP4; TNRC1; Zinc finger protein 384; CAG repeat protein 1; CAS-interacting zinc finger protein; Nuclear matrix transcription factor 4; Nuclear matrix protein 4; Trinucleotide repeat-containing gene 1 protein
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF384

Image



Western Blot analysis of VEC cells using ZNF384 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.