




COL5A2 Antibody

Product Code	CSB-PA006078
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05997
Immunogen	Synthesized peptide derived from the N-terminal region of Human COL5A2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	COL5A2; Collagen alpha-2(V) chain
Immunogen Species	Homo sapiens (Human)
Target Names	COL5A2

Image



Western Blot analysis of HepG2 cells using COL5A2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.