




RSAD2 Antibody

Product Code	CSB-PA006081
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXG1
Immunogen	Synthesized peptide derived from the N-terminal region of Human RSAD2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RSAD2; CIG5; Radical S-adenosyl methionine domain-containing protein 2; Cytomegalovirus-induced gene 5 protein; Viperin; Virus inhibitory protein, endoplasmic reticulum-associated, interferon-inducible
Immunogen Species	Homo sapiens (Human)
Target Names	RSAD2

Image



Western Blot analysis of HepG2 cells using RSAD2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.