




C4A/C4B Antibody

Product Code	CSB-PA006104
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0L4/P0C0L5
Immunogen	Synthesized peptide derived from the Internal region of Human C4a/b.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	C4A; CO4; CPAMD2; Complement C4-A; Acidic complement C4; C3 and PZP-like alpha-2-macroglobulin domain-containing protein 2; C4B; CO4; CPAMD3; Complement C4-B; Basic complement C4; C3 and PZP-like alpha-2-macroglobulin domain-containing protein 3
Immunogen Species	Homo sapiens (Human)
Target Names	C4A/C4B

Image



Western Blot analysis of AD293 cells using C4a/b Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.