




MYOD1 Antibody

Product Code	CSB-PA006107
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15172
Immunogen	Synthesized peptide derived from the Internal region of Human MyoD.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MYOD1; BHLHC1; MYF3; MYOD; Myoblast determination protein 1; Class C basic helix-loop-helix protein 1; bHLHc1; Myogenic factor 3; Myf-3
Immunogen Species	Homo sapiens (Human)
Target Names	MYOD1

Image



Western Blot analysis of NIH-3T3 HeLa cells using MyoD Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.