




IL6R Antibody

Product Code	CSB-PA006108
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08887
Immunogen	Synthesized peptide derived from the Internal region of Human IL-6R α .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL6R; Interleukin-6 receptor subunit alpha; IL-6 receptor subunit alpha; IL-6R subunit alpha; IL-6R-alpha; IL-6RA; IL-6R 1; Membrane glycoprotein 80; gp80; CD126
Immunogen Species	Homo sapiens (Human)
Target Names	IL6R

Image



Western Blot analysis of AD293, 4T1 cells using IL-6R α Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.