




APOC3 Antibody

Product Code	CSB-PA006112
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02656
Immunogen	Synthesized peptide derived from the C-terminal region of Human ApoC-III.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APOC3; Apolipoprotein C-III; Apo-CIII; ApoC-III; Apolipoprotein C3
Immunogen Species	Homo sapiens (Human)
Target Names	APOC3

Image



Western Blot analysis of 293T-UV cells using ApoC-III Polyclonal Antibody