




ALPL Antibody

Product Code	CSB-PA006129
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05186
Immunogen	Synthesized peptide derived from the N-terminal region of Human TNAP.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ALPL; Alkaline phosphatase, tissue-nonspecific isozyme; AP-TNAP; TNSALP; Alkaline phosphatase liver/bone/kidney isozyme
Immunogen Species	Homo sapiens (Human)
Target Names	ALPL

Image



Western Blot analysis of MCF7 cells using TNAP Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.