




GHR Antibody

Product Code	CSB-PA006152
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10912
Immunogen	Synthesized peptide derived from the N-terminal region of Human GHR.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GHR; Growth hormone receptor; GH receptor; Somatotropin receptor
Immunogen Species	Homo sapiens (Human)
Target Names	GHR

Image



Western Blot analysis of SKOV3 cells using GHR Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.