




VIPR1 Antibody

Product Code	CSB-PA006153
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32241
Immunogen	Synthesized peptide derived from the Internal region of Human VPAC1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VIPR1; Vasoactive intestinal polypeptide receptor 1; VIP-R-1; Pituitary adenylate cyclase-activating polypeptide type II receptor; PACAP type II receptor; PACAP-R-2; PACAP-R2; VPAC1
Immunogen Species	Homo sapiens (Human)
Target Names	VIPR1

Image



Western Blot analysis of HepG2 cells using VIPR1 Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.