





CD33 Antibody

Product Code	CSB-PA006175
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20138
Immunogen	Synthesized peptide derived from the Internal region of Human CD33.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD33; SIGLEC3; Myeloid cell surface antigen CD33; Sialic acid-binding Ig-like lectin 3; Siglec-3; gp67; CD33
Immunogen Species	Homo sapiens (Human)
Target Names	CD33

Image



Western Blot analysis of 3T3 cells using CD33 Polyclonal Antibody



Western Blot analysis of NIH-3T3 KB cells using CD33 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.