





CEACAM3/CEACAM6 Antibody

Product Code	CSB-PA006209
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40198/P40199
Immunogen	Synthesized peptide derived from the Internal region of Human CD66c/d.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CEACAM3; CD66D; CGM1; Carcinoembryonic antigen-related cell adhesion molecule 3; Carcinoembryonic antigen CGM1; CD66d; CEACAM6; NCA; Carcinoembryonic antigen-related cell adhesion molecule 6; Non-specific crossreacting antigen; Normal cross-reacting antig
Immunogen Species	Homo sapiens (Human)
Target Names	CEACAM3/CEACAM6

Image



Western Blot analysis of HEPG2 cells using CD66c/d Polyclonal Antibody



Western Blot analysis of HepG2 mouse brain cells using CD66c/d Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.