




SLC3A2 Antibody

Product Code	CSB-PA006233
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08195
Immunogen	Synthesized peptide derived from the C-terminal region of Human CD98.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SLC3A2; MDU1; 4F2 cell-surface antigen heavy chain; 4F2hc; 4F2 heavy chain antigen; Lymphocyte activation antigen 4F2 large subunit; CD98
Immunogen Species	Homo sapiens (Human)
Target Names	SLC3A2

Image



Western Blot analysis of HEB cells using CD98 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.