





PVRL2 Antibody

Product Code	CSB-PA006259
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92692
Immunogen	Synthesized peptide derived from the Internal region of Human Nectin 2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PVRL2; HVEB; PRR2; Poliovirus receptor-related protein 2; Herpes virus entry mediator B; Herpesvirus entry mediator B; HveB; Nectin-2; CD112
Immunogen Species	Homo sapiens (Human)
Target Names	NECTIN2

Image



Western Blot analysis of MOUSE-BRAIN cells using Nectin 2 Polyclonal Antibody



Western Blot analysis of mouse brain cells using Nectin 2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.