




MST1R Antibody

Product Code	CSB-PA006266
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04912
Immunogen	Synthesized peptide derived from the N-terminal region of Human Ron α .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MST1R; PTK8; RON; Macrophage-stimulating protein receptor; MSP receptor; CDw136; Protein-tyrosine kinase 8; p185-Ron; CD136
Immunogen Species	Homo sapiens (Human)
Target Names	MST1R

Image



Western Blot analysis of HeLa cells using Ron α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.