




IL12RB1 Antibody

Product Code	CSB-PA006274
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42701
Immunogen	Synthesized peptide derived from the Internal region of Human IL-12Rβ1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL12RB1; IL12R; IL12RB; Interleukin-12 receptor subunit beta-1; IL-12 receptor subunit beta-1; IL-12R subunit beta-1; IL-12R-beta-1; IL-12RB1; IL-12 receptor beta component; CD212
Immunogen Species	Homo sapiens (Human)
Target Names	IL12RB1

Image



Western Blot analysis of HeLa cells using IL-12Rβ1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.