





CD244 Antibody

Product Code	CSB-PA006275
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZW8
Immunogen	Synthesized peptide derived from the Internal region of Human CD244.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD244; 2B4; Natural killer cell receptor 2B4; NK cell activation-inducing ligand; NAIL; NK cell type I receptor protein 2B4; NKR2B4; h2B4; CD244
Immunogen Species	Homo sapiens (Human)
Target Names	CD244

Image



Western Blot analysis of KB cells using CD244 Polyclonal Antibody



Western Blot analysis of KB cells using CD244 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.