





TNFSF4 Antibody

Product Code	CSB-PA006280
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23510
Immunogen	Synthesized peptide derived from the Internal region of Human Ox40L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFSF4; TXGP1; Tumor necrosis factor ligand superfamily member 4; Glycoprotein Gp34; OX40 ligand; OX40L; TAX transcriptionally-activated glycoprotein 1; CD252
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF4

Image



Western Blot analysis of KB cells using Ox40L Polyclonal Antibody



Western Blot analysis of KB cells using Ox40L Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.