





TNFSF4 Antibody

Product Code	CSB-PA006280
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23510
Immunogen	Synthesized peptide derived from the Internal region of Human Ox40L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFSF4; TXGP1; Tumor necrosis factor ligand superfamily member 4; Glycoprotein Gp34; OX40 ligand; OX40L; TAX transcriptionally-activated glycoprotein 1; CD252
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF4

Image



Western Blot analysis of KB cells using Ox40L Polyclonal Antibody



Western Blot analysis of KB cells using Ox40L Polyclonal Antibody