





NCR3 Antibody

Product Code	CSB-PA006291
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14931
Immunogen	Synthesized peptide derived from the Extracellular and Helical region of Human NKp30.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NCR3

Image



Western Blot analysis of KB cells using NKp30 Polyclonal Antibody



Western Blot analysis of KB cells using NKp30 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.