





LAIR1 Antibody

Product Code	CSB-PA006294
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GTX8
Immunogen	Synthesized peptide derived from the Internal region of Human CD305.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LAIR1; CD305; Leukocyte-associated immunoglobulin-like receptor 1; LAIR-1; hLAIR1; CD305
Immunogen Species	Homo sapiens (Human)
Target Names	LAIR1

Image



Western Blot analysis of KB cells using CD305 Polyclonal Antibody



Western Blot analysis of KB NIH-3T3 cells using CD305 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.