




AREG Antibody

Product Code	CSB-PA006314
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15514
Immunogen	Synthesized peptide derived from the Internal region of Human Amphiregulin.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AREG; SDGF; AREGB; Amphiregulin; AR; Colorectum cell-derived growth factor; CRDGF
Immunogen Species	Homo sapiens (Human)
Target Names	AREG

Image



Western Blot analysis of HBE cells using Amphiregulin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.